

Save running cost and resources

Observe indicator through visual window

No harm to instruments


Easy to pack

Guarantee the sterilization

Durable and lighter



Rapid Sterilization Check
No need to open the container for checking the result of sterilization. You can observe the indicator through the visual window.




Instruction Sticker
Manual record of the sterilization date, expiry date, user and instruments in the container.



Durable Construction
Fix the instruments with silicone bracket in the container.



Locking System
Effective locking and nice feeling while using. Prevent the container opening while transporting and ensure sterilized condition.




Guarantee the Sterilization
We adopt American tyvek 1059B as Container Disk Membrane. Locking system prevents the container opening while transporting and ensures the sterilized condition.

Model	Length (Outer/Inner)	Width (Outer/Inner)	Height (Outer/Inner)	Weight
603013	641/587mm	312/256mm	100/79mm	3.2kg
482010	549/494mm	280/224mm	110/89mm	2.6kg
301010	672/617mm	202/147mm	90/65mm	2.4kg

Standard Accessories




Chemical Indicator Strip



Disk Membrane



Instruction Sticker




Plastic Seals

- ▶ Three kinds of Chemical Indicator Strip & Disk Membrane are available. But, only one kind(Steam/VH₂O₂/ETO) will be supplied as standard accessories.

Accessories Information

For Steam



before use



after use




before use




after use

For VH₂O₂



before use



after use




before use




after use

For ETO



before use



after use



before use



after use

Quality Certificate



ISO 11607



EN 868

Test Report from CDC, China



CASP[®]

Container for Steam/VH₂O₂/ETO
Model: 603013, 482010, 301010



Manufacturer: MuDanJiang Plasma Physical Application Technology Co., Ltd.

Address: JiangNan New Area, Dong'An District, MuDanJiang, HeiLongJiang Province, China

Tel: 82-2-830-4182 Fax: 82-2-830-4183 E-mail: bjderun@naver.com

CASP is our company's registered trademark and it is strictly prohibited to use without our permission www.caspmedical.com